


LRM302212 Armrest Assembly



Item	Part Number	Description	Qty.
1	LRM302212-01	Armrest Pad	2
2	LRM302212-02	Armrest Round Tube	2
3	LRM302212-03	Hex Socket Cap Screws	2
4	LRM302212-04	Washer	2
5	LRM302212-05	Hex Socket Cap Screws	2
6	LRM302212-06	Hex Socket Cap Screws	2
7	LRM302212-07	Armrest Elbow Bracket	2
8	LRM302212-08	Square Cap	2
9	LRM302212-09	Hex Socket Cap Screws	2
10	LRM302212-10	Washer	2
11	LRM302212-11	Spring Washer	2
12	LRM302212-12	Cap Nut	2