


LRM402212 Armrest Assembly



Item	Part Number	Description	Qty.
1	LRM402212-01	Armrest Pad	2
2	LRM402212-02	Armrest Round Tube	2
3	LRM402212-03	Hex Socket Cap Screws	2
4	LRM402212-04	Washer	2
5	LRM402212-05	Hex Socket Cap Screws	2
6	LRM402212-06	Hex Socket Cap Screws	2
7	LRM402212-07	Armrest Elbow Bracket	2
8	LRM402212-08	Square Cap	2
9	LRM402212-09	Hex Socket Cap Screws	2
10	LRM402212-10	Washer	2
11	LRM402212-11	Spring Washer	2
12	LRM402212-12	Cap Nut	2